



Two Tiered Sterilization Container - 25 x 12 x 3cm

This plastic sterilization container has two tiers of silicone slots that will accommodate up to 12 different instruments. Hinged top and bottom lock tight and are perforated. Suitable for autoclave or ethylene oxide sterilization.



This product is discontinued and has limited availability.

General Info

Product Number	20820-00
Category	Instrument Care > Sterilization
EAN	4251909507558



Find 20820-00 on
finescience.com

Technical Info

Shipping Weight	250 g	Net Weight	248 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Plastic		
Autoclavesafe	Yes		
Dimensionstext	25 x 12 x 3 cm		