



Colored Plastic Instrument Case - 17 x 8 x 2cm/Clear



Good for instruments up to 16cm long, this plastic case is safe, secure, and slotted for placement of individual surgical instruments aiding in the prevention of tip damage. Case is autoclavable but should not be used for sterilizing instruments. Instruments not included.

General Info

Product Number	20830-02
Category	Instrument Care > Storage Most Popular > Instrument Care
EAN	4251909507589



Find 20830-02 on
finescience.com

Technical Info

Shipping Weight	70 g	Net Weight	69 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Plastic		
Autoclavesafe	Yes		
Color	Clear		
Dimensionstext	17 x 8 x 2 cm		