



Colored Plastic Instrument Case - 17 x 8 x 2cm/Yellow



Good for instruments up to 16cm long, this plastic case is safe, secure, and slotted for placement of individual surgical instruments aiding in the prevention of tip damage. The case should not be used for sterilizing surgical instruments.

General Info

Product Number	20830-03
Category	Instrument Care > Storage Most Popular > Instrument Care
EAN	4251909507596



Find 20830-03 on
finescience.com

Technical Info

Shipping Weight	70 g	Net Weight	69 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Plastic		
Autoclavesafe	Yes		
Color	Yellow		
Dimensionstext	17 x 8 x 2 cm		