



This sterilization container is made from Ultem polymer to ensure strength and long lasting durability. The transparent perforated lid provides easy identification of your surgical instruments. Silicone mat is included for a safe and secure hold. Trays are also stackable for convenience. Suitable for autoclave or ethylene oxide sterilization.

General Info

Product Number	20840-01
Category	Instrument Care > Sterilization
EAN	4251909507633



Find 20840-01 on
finescience.com

Technical Info

Shipping Weight	145 g	Net Weight	141 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Plastic		
Autoclavesafe	Yes		
Dimensionstext	19 x 7 x 2 cm		