




Plastic Sterilization Case - Double Tier/25 x 15 x 4cm



This sterilization container is made from Ultem polymer to ensure strength and long lasting durability. The transparent perforated lid provides easy identification of your surgical instruments. Silicone mat is included for a safe and secure hold. Trays are also stackable for convenience. Suitable for autoclave or ethylene oxide sterilization.

General Info

Product Number	20842-03
Category	Instrument Care > Sterilization Most Popular > Instrument Care
EAN	4251909507664

Find 20842-03 on
finescience.com

Technical Info

Shipping Weight	805 g	Net Weight	800 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Plastic		
Autoclavesafe	Yes		
Dimensionstext	25 x 15 x 4 cm		