



Stainless Steel Sterilization Container - 30 x 14 x 4cm

This container assists in the sterilization and storage of your fine surgical instruments. It comes with a custom fit silicon mat for keeping the instruments safe during sterilization and storage. Required long-term filters (No. 20850-10) are sold separately. Suitable for autoclave or ethylene oxide sterilization.



General Info

Product Number	20850-00
Category	Instrument Care > Sterilization
EAN	4251909507671

Find 20850-00 on
finescience.com

Technical Info

Shipping Weight	920 g	Net Weight	918 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Aluminum tub with surgical stainless steel (Inox) lid		
Antimagnetic	Yes		
Autoclavesafe	Yes		
Dimensionstext	30 x 14 x 4 cm		