



Filters for Stainless Steel Sterilization Container (Pack of 5)



Use with our Dual Clasping Stainless Steel Sterilization Container No. 20850-00. Filters keep instruments sterile while they remain in the container.

General Info

Product Number	20850-10
Category	Instrument Care > Sterilization
EAN	4251909507688



Find 20850-10 on
finescience.com

Technical Info

Shipping Weight	60 g	Net Weight	55 g
Package Size	5 Pc.	Unit of Measure	PA
Minimum Order Quantity	1		
Alloymaterial	Textile		
Autoclavesafe	Yes		
Dimensionstext	17 x 8 x 2 cm		
Feature	Pack of 5		