



Marking Pen - Fine/Yellow



The Marking Pen was developed for easy, long lasting identification of laboratory animals. The markers contain no harmful chemicals and instead nontoxic pigment approved by the Japanese Ministry of Health for use in food, drugs and cosmetics. It is intended for use on the fur rather than the skin and will last for 6-12 weeks.

General Info

Product Number	24202-03
Category	Lab Accessories > Marking & Identification
EAN	4251909512460



Find 24202-03 on
finescience.com

Technical Info

Shipping Weight	21 g	Net Weight	21 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Autoclavesafe	Yes		
Capacity	10 ml		
Color	Yellow		
Tipdimensions	4.5 x 7.5 mm		
Tips	Fine Tip		