



Marking Pen - Standard/Blue



Capacity: 15ml Tip Size: 15mm x 10mm x 10mm Nontoxic and devoid of any hazardous materials. A single marker can be used to mark 300-500 rats.(depending upon animal size and/or applied area) Once applied, colors remain for 6-12 weeks. Broad and Fine tip Ma

General Info

|                |                                            |
|----------------|--------------------------------------------|
| Product Number | 24202-12                                   |
| Category       | Lab Accessories > Marking & Identification |
| EAN            | 4251909512514                              |



Find 24202-12 on  
[finescience.com](https://finescience.com)

Technical Info

|                        |                 |                 |      |
|------------------------|-----------------|-----------------|------|
| Shipping Weight        | 35 g            | Net Weight      | 35 g |
| Package Size           | 1 Pc.           | Unit of Measure | EA   |
| Minimum Order Quantity | 1               |                 |      |
| Capacity               | 15 ml           |                 |      |
| Color                  | Blue            |                 |      |
| Tipdimensions          | 15 x 10 x 10 mm |                 |      |
| Tips                   | Standard Tip    |                 |      |