



Marking Pen - Standard/Yellow



Capacity: 15ml Tip Size: 15mm x 10mm x 10mm Nontoxic and devoid of any hazardous materials. A single marker can be used to mark 300-500 rats.(depending upon animal size and/or applied area) Once applied, colors remain for 6-12 weeks. Broad and Fine tip Ma

General Info

Product Number	24202-13
Category	Lab Accessories > Marking & Identification
EAN	4251909512521



Find 24202-13 on
finescience.com

Technical Info

Shipping Weight	35 g	Net Weight	35 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Capacity	15 ml		
Color	Yellow		
Tipdimensions	15 x 10 x 10 mm		
Tips	Standard Tip		