



Finger Loop Ear Punch - 1mm Punch Diameter



These ear punches provide users with stability and control when using the numbering system.

General Info

Product Number	24212-01
Category	Rongeurs, Punches & Sampling > Punches
EAN	4251909507749



Find 24212-01 on finescience.com

Technical Info

Shipping Weight	35 g	Net Weight	32 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Autoclavesafe	Yes		
Length	10.5 cm		
Punchdiameter	1 mm		