



Finger Loop Ear Punch - 2mm Punch Diameter



These ear punches provide users with stability and control when using the numbering system. The 1mm ear punch (No. 24212-01) offers greater numbering variations and/or patterns.

General Info

Product Number	24212-02
Category	Rongeurs, Punches & Sampling > Punches
EAN	4251909507756



Find 24212-02 on
finescience.com

Technical Info

Shipping Weight	35 g	Net Weight	30 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Autoclavesafe	Yes		
Length	10.5 cm		
Punchdiameter	2 mm		