



GenoPunch - Stainless Steel & Anodized Aluminum/2mm Punch Diameter



The interchangeable GenoPunch aids in the fast and simple collection of the 2 mm tissue sample for DNA analysis and/or genotyping. GenoPunch includes: (1) GenoPunch (1) 0.5ml Tube Holder (1) 1.6ml Tube Holder (1) Screwdriver (4) Screws FST does not supply microtubes.

General Info

Product Number	24215-02
Category	Rongeurs, Punches & Sampling > Punches Most Popular > Rongeurs, Punches & Sampling
EAN	4251909507763



Find 24215-02 on [finescience.com](https://www.finescience.com)

Technical Info

Shipping Weight	120 g	Net Weight	90 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox) and anodized aluminum		
Autoclavesafe	Yes		
Length	11.5 cm		
Punchdiameter	2 mm		