



Self-Piercing Ear Tags (Black inked) - #1 to #100



Ideal for use in small animals such as mice and rats. These silver colored tags with black numbers are easy to read and are safe for use in salt water. Once applied, the tags cannot be removed. *Applicator sold separately Custom alphanumeric sequences and numbers above 1000 are available by special order. Please call or email our office for more information.

This product is discontinued and has limited availability.

General Info

Product Number	24220-10
Category	Lab Accessories > Marking & Identification
EAN	4251909507787



Find 24220-10 on
[finescience.com](https://www.finescience.com)

Technical Info

Shipping Weight	25 g	Net Weight	21 g
Package Size	100 Pc.	Unit of Measure	PA
Minimum Order Quantity	1		
Alloymaterial	Monel		
Length	0.8 cm		