



Self-Piercing Ear Tags (Laser marked) - #1 to #100



Ideal for use in small animals such as mice and rats, these ear tags come numbered from 001-100. The silver colored ear tags are made of monel and the laser etched numbers provide high visibility making them easy to read. Once applied, tags seal and cannot be removed, providing an easy to use, permanent identification method. These ear tags are safe for use in salt water.

General Info

Product Number	24221-10
Category	Lab Accessories > Marking & Identification Most Popular > Lab Accessories
EAN	4251909511777



Find 24221-10 on
[finescience.com](https://www.finescience.com)

Technical Info

Shipping Weight	25 g	Net Weight	21 g
Package Size	100 Pc.	Unit of Measure	PA
Minimum Order Quantity	1		
Alloymaterial	Monel		
Length	0.8 cm		
Feature	No. 1-100		