



Ear Punch - Stainless Steel/1mm Punch Diameter



The Ear Puncher relieves stress on the wrist during ID marking as it requires only a very gentle pressure to cut a hole. The stainless steel puncher has a sloped edge that provides a cut rather than a punch, which minimizes animal discomfort. The ear puncher is autoclavable but should be lubricated to improve longevity.

General Info

Product Number	24916-01
Category	Rongeurs, Punches & Sampling > Punches
EAN	4251909512118



Find 24916-01 on
finescience.com

Technical Info

Shipping Weight	75 g	Net Weight	74 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Autoclavesafe	Yes		
Length	14 cm		
Punchdiameter	1 mm		