



Insect Pins - Black Enameled/Size 2



With nylon heads.

General Info

Product Number	26000-45
Category	Pins, Probes & Hooks > Pins
EAN	4251909507855



Find 26000-45 on
finescience.com

Technical Info

Shipping Weight	10 g	Net Weight	8 g
Package Size	100 Pc.	Unit of Measure	PA
Minimum Order Quantity	1		
Alloymaterial	Black Enameled Steel		
Length	4 cm		
Roddiameter	0.45 mm		
Size	2		
Tipdiameter	0.03 mm		