



Insect Pins - Stainless Steel/Size 6



With nylon heads.

General Info

Product Number	26001-65
Category	Pins, Probes & Hooks > Pins
EAN	4251909507992



Find 26001-65 on
finescience.com

Technical Info

Shipping Weight	15 g	Net Weight	8 g
Package Size	100 Pc.	Unit of Measure	PA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	4 cm		
Roddiameter	0.65 mm		
Size	6		
Tipdiameter	0.03 mm		