



Insect Pins - Stainless Steel/Size 7



With nylon heads.

General Info

Product Number	26001-70
Category	Pins, Probes & Hooks > Pins
EAN	4251909508005



Find 26001-70 on
finescience.com

Technical Info

Shipping Weight	20 g	Net Weight	8 g
Package Size	100 Pc.	Unit of Measure	PA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	5.5 cm		
Roddiameter	0.7 mm		
Size	7		
Tipdiameter	0.03 mm		