



Minutien Pins - Stainless Steel/0.1mm Diameter/0.0125mm Tip



General Info

Product Number	26002-10
Category	Pins, Probes & Hooks > Pins Most Popular > Pins, Probes & Hooks
EAN	4251909508012



Find 26002-10 on
finescience.com

Technical Info

Shipping Weight	10 g	Net Weight	8 g
Package Size	500 Pc.	Unit of Measure	PA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	1.2 cm		
Roddiameter	0.1 mm		
Tipdiameter	0.0125 mm		