



Moria Nickel Plated Pin Holder - 12cm



Works with all pins and tungsten needles sold by FST.

General Info

Product Number	26016-12
Category	Pins, Probes & Hooks > Holders Most Popular > Pins, Probes & Hooks
EAN	4251909508098



Find 26016-12 on
finescience.com

Technical Info

Shipping Weight	15 g	Net Weight	12 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Nickel-plated		
Length	12 cm		
Jawopeningdiameter	0 to 1 mm		