



Nickel Plated Pin Holder - 17cm



Works with all pins and tungsten needles sold by FST.

General Info

Product Number	26018-17
Category	Pins, Probes & Hooks > Holders Most Popular > Pins, Probes & Hooks
EAN	4251909508104



Find 26018-17 on
finescience.com

Technical Info

Shipping Weight	15 g	Net Weight	11 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Aluminum handle with nickel-plated chuck		
Antimagnetic	Yes		
Length	17 cm		
Jawopeningdiameter	0 to 1 mm		