



Live Insect Forceps - Sharp/Straight/10cm



General Info

Product Number	26029-10
Category	Forceps & Tweezers > Smooth
EAN	4251909508111



Find 26029-10 on
finescience.com

Technical Info

Shipping Weight	5 g	Net Weight	1 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	10 cm		
Tipdimensions	0.9 x 0.15 mm		
Tips	Narrow		
Tipshape	Straight		