



Three activities can be timed up or down simultaneously or individually. Time can be displayed either in 12 or 24 hours and features a memory function. Maximum countdown time is 19 hours, 59 minutes, 59 seconds. Minimum time interval is 1 second. The timer housing includes a powerful magnet for mounting.

General Info

Product Number	27000-01
Category	Lab Accessories
EAN	4009816014757



Find 27000-01 on
finescience.com

Technical Info

Shipping Weight	115 g	Net Weight	110 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Dimensionstext	7 x 8.5 x 1.5 cm		
Lcddisplay	4.5 x 4.5 cm		
Voltagepowersource	AA battery (not included)		