



Self-Sealing Sterilization Pouches - 10 x 17.6cm



Preformed pouches with self adhering adhesive tape are blue tinted on one side for visualization of contents, and paper-backed on the other side. Pouches will work with either steam autoclave or ethylene oxide sterilization.

- Advantages:
- Protects instruments during sterilization.
 - Provides sterile storage until use.
 - Write on pouches for easier identification during sterilization.

This product is discontinued and has limited availability.

General Info

Product Number	28000-07
Category	Instrument Care > Sterilization
EAN	4251909508142



Find 28000-07 on
finescience.com

Technical Info

Shipping Weight	285 g	Net Weight	2.8 g
Package Size	100 Pc.	Unit of Measure	PA
Minimum Order Quantity	1		
Alloymaterial	Medical grade paper		
Length	21 cm		
Width	10.5 cm		