



Self-Sealing Sterilization Pouches - 11.5 x 22.5 cm



Preformed pouches with self adhering adhesive tape are blue tinted on one side for visualization of contents, and paper-backed on the other side. Pouches will work with either steam autoclave or ethylene oxide sterilization.

- Advantages:
- Protects instruments during sterilization.
  - Provides sterile storage until use.
  - Write on pouches for easier identification during sterilization.

General Info

|                |                                 |
|----------------|---------------------------------|
| Product Number | 28000-11                        |
| Category       | Instrument Care > Sterilization |
| EAN            | 4251909508159                   |



Find 28000-11 on  
[finescience.com](https://finescience.com)

Technical Info

|                        |                                            |                 |       |
|------------------------|--------------------------------------------|-----------------|-------|
| Shipping Weight        | 458 g                                      | Net Weight      | 4.5 g |
| Package Size           | 100 Pc.                                    | Unit of Measure | PA    |
| Minimum Order Quantity | 1                                          |                 |       |
| Alloymaterial          | Medical grade paper                        |                 |       |
| Autoclavesafe          | Yes                                        |                 |       |
| Length                 | 22.5 cm                                    |                 |       |
| Width                  | 11.5 cm                                    |                 |       |
| Dimensionstext         | Outer: 13.5 x 28 cm  Inner: 11.5 x 22.5 cm |                 |       |