



Self-Sealing Sterilization Pouches - 17 x 30 cm



Preformed pouches with self adhering adhesive tape are blue tinted on one side for visualization of contents, and paper-backed on the other side. Pouches will work with either steam autoclave or ethylene oxide sterilization.

- Advantages:
- Protects instruments during sterilization.
 - Provides sterile storage until use.
 - Write on pouches for easier identification during sterilization.

General Info

Product Number	28000-14
Category	Instrument Care > Sterilization
EAN	4251909508166



Find 28000-14 on
finescience.com

Technical Info

Shipping Weight	828 g	Net Weight	8.2 g
Package Size	100 Pc.	Unit of Measure	PA
Minimum Order Quantity	1		
Alloymaterial	Medical grade paper		
Autoclavesafe	Yes		
Length	30 cm		
Width	17 cm		
Dimensionstext	Outer: 19 x 36 cm Inner: 17 x 30 cm		