



Autoclave Tape - Sterilized/2.6cm



Tape will permanently indicate the word STERILIZED or diagonal stripes after exposure to steam autoclave temperature between 121C and 130C during a 40-60 minute time period. Use with sterilization pouches, sterilization wraps and sterilization cases.

This product is discontinued and has limited availability.

General Info

Product Number	28001-01
Category	Instrument Care > Sterilization
EAN	4251909509798



Find 28001-01 on
finescience.com

Technical Info

Shipping Weight	212 g	Net Weight	210 g
Package Size	1 Pc.	Unit of Measure	RO
Minimum Order Quantity	1		
Autoclavesafe	Yes		
Lengthofroll	550 cm		
Tapewidth	2.6 cm		