



Autoclave Tape - Striped/2cm



Tape will permanently indicate the word STERILIZED or diagonal stripes after exposure to steam autoclave temperature between 121C and 130C during a 40-60 minute time period. Use with sterilization pouches, sterilization wraps and sterilization cases.

General Info

Product Number	28001-02
Category	Instrument Care > Sterilization
EAN	4251909508180



Find 28001-02 on
finescience.com

Technical Info

Shipping Weight	130 g	Net Weight	128 g
Package Size	1 Pc.	Unit of Measure	RO
Minimum Order Quantity	1		
Dimensionstext	1,90 cm x 55 m		
Lengthofroll	55 m		