



Lens for Adjustable Gooseneck Stand - 2.8x



Use with Adjustable Gooseneck Stand (No. 28010-00 - Discontinued).

This product is discontinued and has limited availability.

General Info

Product Number	28010-28
Category	Lab Accessories > Magnifiers
EAN	4251909509828



Find 28010-28 on
[finescience.com](https://www.finescience.com)

Technical Info

Shipping Weight	200 g	Net Weight	185 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Lenssize	100 x 75mm		
Lensshape	Aspheric		
Magnification	2.8x		