



Eschenbach Top-Lite LED

Six LEDs are encased along the black rod for high contrast illumination. Three universal clamps are included for mounting. Use with Eschenbach Tilttable Magnifier Stand (No. 28030-38: Discontinued) and 2x Magnification Lens for Adjustable Gooseneck Stand (No. 28010-02).



This product is discontinued and has limited availability.

General Info

Product Number	28015-00
Category	Lab Accessories > Magnifiers
EAN	4251909509866



Find 28015-00 on
[finescience.com](https://www.finescience.com)

Technical Info

Shipping Weight	300 g	Net Weight	278 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Cablelength	2 m		