



Disposable Sterilization Wraps - 38 x 38cm

These biodegradable cellulose fiber wraps are resistant to tears and conform easily to maintain a breathable yet impervious barrier to airborne bacteria, saline, blood and water. Wrapping your instrument case securely will ensure the sterility of your instruments until ready for use. Available in two sizes.



General Info

Product Number	28015-15
Category	Instrument Care > Sterilization
EAN	4251909508197



Find 28015-15 on
finescience.com

Technical Info

Shipping Weight	335 g	Net Weight	332 g
Package Size	25 Pc.	Unit of Measure	PA
Minimum Order Quantity	1		
Alloymaterial	Biodegradable Cellulose Fiber		
Autoclavesafe	Yes		
Dimensionstext	38 x 38 cm		