



Eschenbach varioLEDflex Magnifier - 2.5x



The Vario LED Flex ensures ease and precision with its smooth flexible neck and tiltable magnifying lens. This magnifiers innovative features provide a spacious working platform suitable for any application. Mounting clamp provides uncompromised workspace due to its sturdy slim design. The magnification lens is illuminated with shadow-free full or half ring optimum bright LEDs.

General Info

Product Number	28015-25
Category	Lab Accessories > Magnifiers
EAN	4251909508203



Find 28015-25 on
[finescience.com](https://www.finescience.com)

Technical Info

Shipping Weight	1305 g	Net Weight	1170 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Diopter	6D		
Illumination	SMD LED		
Lenssize	76 mm		
Lensshape	Aspheric		
Lensmaterial	PXM® lightweight lens with patented Ceratec®scratch resistant coating		
Magnification	2.5x		
Voltagepowersource	120 V & 230 V		
Workingdistance	Variable		