



Headband Magnifier - 2.5x



Precision ground and polished optical glass lenses. The padded headband is lightweight and fully adjustable. Available in 2.5x, 2.75x and 3.5x magnification. Can be worn with spectacles.

General Info

Product Number	28030-02
Category	Lab Accessories > Magnifiers
EAN	4251909508234



Find 28030-02 on
[finescience.com](https://www.finescience.com)

Technical Info

Shipping Weight	280 g	Net Weight	150 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Various		
Magnification	2.5x		
Workingdistance	20 cm (8 inch)		