



Headband Magnifier - 2.75x



Precision ground and polished optical glass lenses. The padded headband is lightweight and fully adjustable. Available in 2.5x, 2.75x and 3.5x magnification. Can be worn with spectacles.

General Info

|                |                              |
|----------------|------------------------------|
| Product Number | 28030-03                     |
| Category       | Lab Accessories > Magnifiers |
| EAN            | 4251909508241                |



Find 28030-03 on  
[finescience.com](https://www.finescience.com)

Technical Info

|                        |                |                 |       |
|------------------------|----------------|-----------------|-------|
| Shipping Weight        | 260 g          | Net Weight      | 150 g |
| Package Size           | 1 Pc.          | Unit of Measure | EA    |
| Minimum Order Quantity | 1              |                 |       |
| Alloymaterial          | Various        |                 |       |
| Magnification          | 2.75x          |                 |       |
| Workingdistance        | 15 cm (6 inch) |                 |       |