



Eschenbach maxDETAIL Spectacles - 2x

These spectacles are portable, lightweight, and easily adjustable. Designed with frosted objective lens housing which eliminates distracting reflections. Protective case is included. Shown with Eschenbach Clip-On Light, sold separately.



General Info

Product Number	28030-20
Category	Lab Accessories > Magnifiers MedTech Line
EAN	4048347515753



Find 28030-20 on
[finescience.com](https://www.finescience.com)

Technical Info

Shipping Weight	230 g	Net Weight	114 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Fieldofview	110 mm (4 inch)		
Magnification	2x		
Workingdistance	40 cm (16 inch)		