



Eschenbach ridoMED Spectacles with Band - 2.5x

These magnifier spectacles allow full adjustments on almost all parts. Headband keeps spectacles steady and in place (sold separately). Leather case included.



General Info

Product Number	28040-30
Category	Lab Accessories > Magnifiers
EAN	4251909508272



Find 28040-30 on
finescience.com

Technical Info

Shipping Weight	140 g	Net Weight	97 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Plastic frame reinforced with fiberglass		
Magnification	2.5x		
Workingdistance	35 cm (14 inch)		