



Galilean Flip-up Spectacles - 3x



These loupes are lightweight and easy to use. Made from the finest optical coated lenses which deliver excellent clarity and high resolution. The magnifier loupes are mounted onto a strong yet flexible wraparound safety frame with a headband for stability and comfort. Ability to flip loupes up for use as safety glasses. Comes with mini screwdriver for fine adjustments, micro cloth, and storage case. Available in 3x and 3.5x magnification.

General Info

Product Number	28050-30
Category	Lab Accessories > Magnifiers
EAN	4251909508289



Find 28050-30 on
[finescience.com](https://www.finescience.com)

Technical Info

Shipping Weight	270 g	Net Weight	262 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Various		
Magnification	3x		
Workingdistance	40 cm (16 inch)		