



Economical and fully interchangeable headband magnifier with LED light. Includes: Multiple Lens Magnifications (1.0x, 1.5x, 2.0x, 2.5x, 3.5x) LED Light Source Tilting Functionality Cushioned Nose Pads Adjustable Headband Interchangeable Lens Frames Requires 3 AAA Batteries (Batteries not included)

### General Info

Product Number	28080-10
Category	Lab Accessories > Magnifiers MedTech Line
EAN	4251909508302



Find 28080-10 on  
[finescience.com](https://www.finescience.com)

### Technical Info

Shipping Weight	315 g	Net Weight	213 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Autoclavesafe	Yes		
Fieldofview	Variable		
Magnification	1x, 1.5x, 2x, 2.5x and 3.5x		
Weeeno	FST		
Workingdistance	Variable		