



Repair Kit for Dumont Forceps



Realign and repair Dumont forceps with this convenient kit. Bent forceps tips can be re-straightened with the chain nose pliers. The sharpening stone has two grooves, 1.6 mm and 0.8 mm, to assist in resharpenering the forceps tips. O-rings clamp the forceps shut to ensure that both shanks are equally aligned and sharpened. Instrument oil assists in producing a smooth sharp edge. Detailed instructions are included.

- Content of Kit:
- 1 Sharpening Stone [29008-22]
 - 1 Bottle of Instrument Oil (50ml) [29055-00]
 - 1 Flat Nose Pliers - not sold seperately
 - 1 Pack of 10 O-Rings [11200-00]

General Info

Product Number	29000-00
Category	Instrument Care > Maintenance Kits Most Popular > Instrument Care
EAN	4251909508326



Find 29000-00 on finescience.com

Technical Info

Shipping Weight	415 g	Net Weight	400 g
Package Size	1 Pc.	Unit of Measure	SET
Minimum Order Quantity	1		
Alloymaterial	Various		