




Sharpening Stone - Square



Natural Arkansas oil stones feature unequalled fineness, hardness and structural density. When used with instrument oil, they will produce the sharpest cutting edge. Item (No. 29008-01) shown here.

General Info

Product Number	29008-01
Category	Instrument Care > Maintenance
EAN	4251909508340



Find 29008-01 on
finescience.com

Technical Info

Shipping Weight	15 g	Net Weight	9 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Arkansas stone		
Dimensionstext	6 mm x 9 cm		
Shape	Square		