



Sharpening Stone - Triangle



Natural Arkansas oil stones feature unequalled fineness, hardness and structural density. When used with instrument oil, they will produce the sharpest cutting edge.

General Info

Product Number	29008-12
Category	Instrument Care > Maintenance
EAN	4251909508357



Find 29008-12 on
finescience.com

Technical Info

Shipping Weight	10 g	Net Weight	4 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Arkansas stone		
Dimensionstext	6 mm x 9 cm		
Shape	Triangle		