



Sharpening Stone for Fine Forceps



This natural Arkansas stone assists in resharpener the tips of Dumont forceps. The stone has two grooves, 0.8mm and 1.6mm, and comes in cedar box with ten O-rings.

General Info

Product Number	29008-22
Category	Instrument Care > Maintenance
EAN	4251909508364



Find 29008-22 on
finescience.com

Technical Info

Shipping Weight	230 g	Net Weight	226 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Arkansas stone		
Dimensionstext	12.5 x 6.5 x 3 cm		
Shape	Flat with two Grooves		