

These are highly accurate stage micrometers for use with transmitted light in both compound and stereo microscopes. Comes with a plastic case.



This product is discontinued and has limited availability.

General Info

Product Number	29025-01
Category	Lab Accessories > Measurement
EAN	4251909508388



Find 29025-01 on
[finescience.com](https://www.finescience.com)

Technical Info

Shipping Weight	20 g	Net Weight	5 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Autoclavesafe	Yes		
Dimensionstext	76 x 26 mm		
Increments	10 mm with 100 divisions (0.1 mm)		
Maximummeasurement	10 mm		