



Stage Micrometer - 100 divisions of 0.01mm

These are highly accurate stage micrometers for use with transmitted light in both compound and stereo microscopes. Comes with a plastic case.



General Info

Product Number	29025-02
Category	Lab Accessories > Measurement
EAN	4251909508395



Find 29025-02 on
finescience.com

Technical Info

Shipping Weight	25 g	Net Weight	5 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Autoclavesafe	Yes		
Dimensionstext	76 x 26 mm		
Increments	1 mm with 100 divisions (0.01 mm)		
Maximummeasurement	1 mm		