



Loupe with Micrometer and LED

High quality loupe features a removable micrometer and ultra-bright LED lights surrounding the circumference for a defined image and clear contrast.



General Info

Product Number	29026-10
Category	Lab Accessories > Magnifiers
EAN	4251909508401



Find 29026-10 on
[finescience.com](https://www.finescience.com)

Technical Info

Shipping Weight	160 g	Net Weight	141 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Dimensionstext	38 x 38 x 40 mm		
Fieldofview	30 mm		
Illumination	Ultrabright LED		
Magnification	10x		
Workingdistance	10 mm (without scale)		