



Silicone Lubricant



This stiff non-melting silicone base maintains its consistency from -40C to 200C (40F to 392F) with less than 2 evaporation after 24 hours of exposure to 150C (302F). Typical applications include lubrication of stopcock rotating parts, surgical instruments and other lab items. Seals vacuum and low pressure systems. Compatible with glass, plastics, metals, synthetic rubber gaskets and seals.

This product is discontinued and has limited availability.

General Info

Product Number	29051-45
Category	Instrument Care > Maintenance
EAN	4251909508425



Find 29051-45 on
finescience.com

Technical Info

Shipping Weight	25 g	Net Weight	20 g
Package Size	15 ml	Unit of Measure	MLT
Minimum Order Quantity	1		
Autoclavesafe	Yes		
Volume	15 ml		