




Instrument Tip Protectors - Assorted Pack



These tip protectors are a convenient, color-coded way to protect surgical instrument tips and edges. They provide exceptional protection and are vented to prevent rusting from condensation. The Assorted Pack (Item No. 29060-00) contains five of each size.

General Info

Product Number	29060-00
Category	Instrument Care > Storage Most Popular > Instrument Care
EAN	4251909508449



Find 29060-00 on
finescience.com

Technical Info

Shipping Weight	25 g	Net Weight	22 g
Package Size	30 Pc.	Unit of Measure	PA
Minimum Order Quantity	1		
Alloymaterial	Plastisol		
Autoclavesafe	Yes		
Color	Assorted		
Dimensionstext	Assorted		