



Instrument Tip Protectors - 3.2 x 25.4mm/Clear



These tip protectors are a convenient, color-coded way to protect surgical instrument tips and edges. They provide exceptional protection and are vented to prevent rusting from condensation. The Assorted Pack (Item No. 29060-00) is 5 of each of the 6 available sizes.

General Info

Product Number	29060-04
Category	Instrument Care > Storage
EAN	4251909508487



Find 29060-04 on
[finescience.com](https://www.finescience.com)

Technical Info

Shipping Weight	10 g	Net Weight	8 g
Package Size	25 Pc.	Unit of Measure	PA
Minimum Order Quantity	1		
Alloymaterial	Plastisol		
Autoclavesafe	Yes		
Color	Clear		
Dimensionstext	3.2 x 25.4 mm		