



Instrument Tip Protectors - 4.8 x 25.4mm/Clear



These tip protectors are a convenient, color-coded way to protect surgical instrument tips and edges. They provide exceptional protection and are vented to prevent rusting from condensation. The Assorted Pack (Item No. 29060-00) is 5 of each of the 6 available sizes.

General Info

Product Number	29060-05
Category	Instrument Care > Storage
EAN	4251909508494



Find 29060-05 on
[finescience.com](https://www.finescience.com)

Technical Info

Shipping Weight	20 g	Net Weight	17 g
Package Size	25 Pc.	Unit of Measure	PA
Minimum Order Quantity	1		
Alloymaterial	Plastisol		
Autoclavesafe	Yes		
Color	Clear		
Dimensionstext	4.8 x 25.4 mm		